


Cut out these templates to select the size that fits best.

Medium fits most women, Large fits most men.

step 1 print at 100% 




step 3
cut along
this line

step 2 cut along this line

 **ZSHIELD** r e v e a l

sizing guide

Note : people with very large or small faces may not be able to achieve a good fit with either a Large or Medium mask.



top edge below bridge of nose
bottom edge above point of chin
does not touch bottom lip



top edge at bridge of nose
bottom edge below point of chin



top edge below bridge of nose
bottom edge above point of chin
does not touch bottom lip



top edge at bridge of nose
touches lower lip